SELF-COLLECTION OF VAGINAL SAMPLE RECOMMENDATIONS

		: ر_
Birthname:		I
Surname:		For use by the
Forename:		laboratory only
Date of birth: / /		
CNS no:	Sampla	ر_ collection
E-mail:@	 Date:	
. 00/+	Time:	
O1. Wash your hands before you begin		02. Open the bag and take out the swab From this point on, do not place the swab on a surface and do not touch the cotton
O3. Hold the swab as shown in the picture and stand or sit on the toilet, with your legs apart		04. Spread your lips and gently introduce the swab into your vagina over 3 to 5 cm Rub the sides of the vagina with the swab for 10 seconds
05. Unscrew the tube cap and insert the swab Do not remove the liquid in the tube!		06. Break the stem of the swab by pressing the pre-formed breaking area (pink) against the wall of the tube
07. Put the cap back on the tube and write your surname and forename Complete the identification box at the top of this form	I The second sec	08. Hand over the sample and form in one of our centres within the storage time indicated below
SAMPLE STORAGE +20°C Max 24 hours		Save time, pre-register
BIONEXT LABORATOIRE D'ANALYSES MÉDICALES F_PREA_126_001 WWW.K	L-336	Je du Château d'Eau 4 Leudelange nbourg ■ (+352) 27 321 (€ (+352) 27 321 201 (€ (+352) 27 321 201 (€ info@bionext.lu