



Forname

Birthname

Surname

Date of birth

CNS N°

@
 E-mail

+
 Phone

For use by the
laboratory only

URINE COLLECTION

Date

Time

Time of passage in tubes

After local cleaning: OUI NON

Collection method:

- 2nd part of stream
 1st part of stream
 Evacuation catheter
 Indwelling catheter
 Bag
 Other

CLINICAL INFORMATION

Please check boxes as appropriate

- | | |
|---|--|
| <input type="checkbox"/> Pain or burning while urinating | <input type="checkbox"/> No clinical signs |
| <input type="checkbox"/> Frequent or uncontrolled need to urinate | <input type="checkbox"/> Check-up before surgical procedure urodynamic exploration |
| <input type="checkbox"/> Fever | <input type="checkbox"/> Kidney transplant |
| <input type="checkbox"/> In chemotherapy | <input type="checkbox"/> Pyelonephritis |
| <input type="checkbox"/> Pregnancy
Estimated due date: _____ | <input type="checkbox"/> Recurrent cystitis |

Did you take antibiotics?

Before collection

- No Yes If yes, which one? _____

Last taken on: /
 DD MM

After collection

- No Yes If yes, which one? _____



URINE COLLECTION RECOMMENDATIONS FOR URINE CULTURE



If possible, collect the first urine of the morning or 4 hours after you last urinated.




Do not touch the integrated tip (inside of the lid).

01. Wash your hands and clean the area around your genitals **with the supplied wipe**.




02. Urinate:
1st part in the toilet
2nd part in the container
Close the container.



03. Lift the white label on the lid and place the **beige** tube with the plug side into the orifice of the container and press down firmly.




04. Repeat with the **khaki** tube. Replace the white label on the hole in the container lid.



05. Turn the tubes over 8 to 10 times.



06. Write your name and **forname, on the 2 tubes and the jar**. Complete the information on the back of this sheet.



07. Insert the tubes, container and the sheet into the bag and seal it.



Finding a BIONEXT centre

08. Take the bag to one of our centres within the storage time indicated below.

URINE CONSERVATION

IN TUBES

+20°C Max 24h

IN CONTAINER

+4°C Max 12h or +20°C Max 2h

Save time
pre-register

